



Cavex is committed to the sales of sophisticated, total systems whereby each individual product is part of a greater whole. Through this so-called systems thinking Cavex provides complete solutions that enable dental surgeons to significantly increase the degree of service to their patients. Your impression is our concern

**The Cavex VacuFormer System:
From impression to the desired splint/mouth
guard/bleaching tray in 1 session..!**

The Cavex VacuFormer System is a complete vacuum thermo forming system comprising the Cavex VacuFormer, deep-drawing plates, a unique plaster composite and accessories.

The Cavex VacuFormer System enables the dentist to optimize the treatment of and service to the patient. Bleaching trays, splints and mouth protectors can be prepared during patient visits without the intervention of the dental laboratory. A detailed protocol, illustrated with photographs, enables anyone to master this technique.





Norsk Orthoform
DENTALARTIKLER



Depot as

www.norskorthoform.no
Telefon 22 76 01 40

The Cavex VacuFormer System: 1 complete system for all indications

● **Cavex
VacuFormer**

The Cavex VacuFormer is a “no-nonsense” all-round thermal vacuum former. This robust machine has an electric, high performance vacuum motor suitable for every conventional dental deep-drawing application. The Cavex VacuFormer is easy to operate and is protected against overheating.

● **Cavex
Instant Stone**

Cavex Instant Stone is a unique plaster composite that can cast a perfect plaster mould in just 10 minutes. Thanks to a sophisticated blending system, Cavex Instant Stone is very clean and easy to use.

● **Cavex Mouth
Protector**

The Cavex mouth protectors are professional 3,8mm dual layer and 5mm triple layer laminates with a combined hard and soft layer. The soft layer has provides very good shock absorption and the hard layer spreads the force of impact and protects the teeth, resulting in a maximum protection of all teeth.

● **Cavex Bruxism**

The Cavex Bruxism laminates are hard, 2mm thick, transparent laminates which can be used to make a retainer preventing bruxism. This retainer can be used in a bruxism treatment plan to prevent further wear of the dentition.

● **Cavex Bleach**

The Cavex Bleach laminates are made of a soft, clear material that is easy to cut and finish. The soft material is comfortable to the gingival area and ideal for making bleaching trays.

● **Cavex Splint**

The Cavex Splint laminates are made of a transparent and rigid material which is ideal for splints, surgical trays, orthodontic retainers, etc. This unique, clear material makes the Cavex Splint invisible when worn in the mouth.

● **Cavex
VacuFormerSystem
Starterkit**

The Starter Kit consists of all equipment and materials that are required for making ready-made deep-drawing applications: 1 Cavex VacuFormer + 18 pots Instant Stone + deep-drawing plates (50 × Splint + 25 × Bleach + 12 × Mouth Protector) + Scissors + 6 Base trays + protocols per indication.

● **Cavex Bleach Tray
Spacer**

Cavex Bleach Tray Spacer is a blue coloured flowable resin especially designed for creating reservoirs in bleach trays.

For further information, please contact Cavex: +31 23 530 77 00 or www.cavex.nl



Assortment:

VF001	Cavex Vacuformer Starterkit	
VT010	Cavex Instant Stone powder liquid	72 sets
VT022	Cavex Bleach/clear/1 mm*	25 sheets
VT030	Cavex Bleach Heavy/clear/1.5 mm*	25 sheets
VT020	Cavex Mouth Protector/red, white & blue/3.8 mm*	12 sheets
VT034	Cavex Mouth Protector/red/3.8 mm*	12 sheets
VT035	Cavex Mouth Protector/blue/3.8 mm*	12 sheets
VT036	Cavex Mouth Protector/clear/3.8 mm*	12 sheets
VT031	Cavex Mouth Protector X-Pro/green/5 mm*	12 sheets
VT032	Cavex Mouth Protector X-Pro/black/5 mm*	12 sheets
VT033	Cavex Mouth Protector X-Pro/clear/5 mm*	12 sheets
VT021	Cavex Splint/clear/0.6 mm*	50 sheets
VT026	Cavex Splint X-heavy/clear/1.5 mm*	25 sheets
VT027	Cavex Bruxism/clear/2 mm*	25 sheets
	* width x length = 12.7 x 12.7 cm	
VT025	Cavex Bleach Tray Spacer	2 x 2 ml syringe + 10 dispensing tips

